## en - Medical Device - Instructions for Use

Intended purpose: To tone, shape, strengthen or rehabilitate muscle groups in fingers, hand and forearms. Intended user: lay person.

Indications: Medical conditions resulting in loss of strength, flexibility and coordination of fingers, hand and forearm muscles. Contra-indications: In case of unusual swelling or discomfort, use should be discontinued.

Precautions:

- Consult a healthcare professional before use;
- Read this manual carefully.

General information:

- Exercises should be done at a slow, comfortable pace;
- Stop using when one or more springs wear out;
- Report any serious incident (death, serious deterioration of health, serious public health threat) to the Manufacturer and your Competent Authority;
- Latex-free.

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# MoVeS

## FLEX-ION

Approved by Certified Hand Therapists

Resistance level: per finger / per hand	extra light	0,35 / 1,1 kg	02-000101	05420063000925
Widerstandskraft: pro finger / gesamtes	light	0,7 / 2,3 kg	02-000102	05420063000932
Niveau de résistance: par doigt / sur l'ensemble	medium	1,4 / 4,5 kg	02-000103	05420063000949
Weerstandsniveau: per vinger / geheel	heavy	2,3 / 7,3 kg	02-000104	05420063000956
Nivel de resistencia: por dedo / sistema completa	extra heavy	3,2 / 10,4 kg	02-000105	05420063000963
Livello di resistenza: per singolo gito / intero apparecchio	super heavy	4,1 / 14,1 kg	02-000106	05420063000970

Daily use conditions and rehabilitates the hands, wrists, and forearms by strengthening and toning the muscles and tendons.

Made in South Korea

UK CA UK REP  
Qarad UK Ltd.  
8 Northumberland Ave  
Westminster, London WC2N 5BY  
United Kingdom

Exercise program enclosed

WWW.MOVES-YOU.COM

MVS In Motion  
Westdijk 150  
2830 Tiselt (Willebroek)  
Belgium -Europe  
www.mvs-in-motion.com  
Find us on &

# MoVeS

## FLEX-ION

Approved by Certified Hand Therapists

EXERCISE EACH FINGER INDIVIDUALLY, OR THE ENTIRE HAND

For additional exercises, developed by a licensed hand therapist, please visit [www.moves-you.com](http://www.moves-you.com)

**en** - Hook Position

Place finger tips on individual buttons and position ergonomic palm bar with the hook over the web space between the thumb and index finger. Keeping finger tips slightly bent, flex all fingers toward the center of the unit.

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**en** - Tip-to- Tip Pinch

Place the tip of the thumb on the center of the ergonomic palm bar. Position the finger to be exercised on one of the center buttons. Keeping both tips slightly bent, flex the thumb and the finger being exercised toward the center of the unit at the same time.

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**en** - Gross Grasp

Place the middle of the fingers on the buttons with the finger tips wrapped over the top. Positioning the hook of the ergonomic palm bar over the web space between the thumb and index finger, flex the thumb and all fingers toward the center of the unit as if making a fist.

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**en** - Fingertip Flexion

Place all finger tips on their own button, with the tip of the thumb on the center of the ergonomic base. Keeping all finger tips bent, flex the thumb and all the fingers toward the center of the unit at the same time.

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**en** - Key Pinch

Place the hook of the ergonomic palm bar on the inside of the bent index finger. With the tip of the thumb slightly bent, flex the thumb on the first button.

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**For additional exercises, developed by a licensed hand therapist, please visit [www.moves-you.com](http://www.moves-you.com)**

**en** - Trigger Pinch

Loosely hold the unit vertically with the buttons toward the fingers and the hook of the ergonomic palm bar positioned over the web space between the thumb and index finger. With the middle of the index finger, press the second button as if squeezing a trigger.

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