# Standard Knock Down Pedal Exerciser

## Standard Knockdown Pedal Exerciser For upper & lower body exercise

### en - Medical Device - Instructions for Use: Lightweight pedal exerciser for arms and legs

Intended purpose: Prevention, reduction, and improvement of the consequences of loss or lack of movement through exercise therapy for arms and legs. Intended user: Lay person.

**Indications:** Prevention, reduction, and improvement of the consequences of loss or lack of movement, mainly in the following medical conditions: neuromuscular, orthopedic, cardiovascular, metabolic and other medical conditions resulting in limitations in movement or becoming confined to bed. **Contra-indications:** Use should be discontinued in case of unusual swelling, discomfort, dizziness, severe muscle or joint pain or pain in the chest.

**Precautions:**

- Consult a healthcare professional before use;
- Use on a flat, stable, non-slip surface;
- Do not wear loose clothing while exercising;
- During leg exercises: always be seated, do not sit in a chair with casters, never stand on the pedal exerciser;
- Read this manual carefully.

**General information:**

- Exercises should be done at a slow, comfortable pace;
- Use the resistance knob to adjust resistance to suit patient's ability;
- Designed as a home use exercise device;
- Do not use when device is damaged.

Report any serious incident (death, serious deterioration of health, serious public health threat) to the Manufacturer and your Competent Authority.

WWW.MOVES-YOU.COM

LOT YYY-MM-DD

MVS IN MOTION

YYY-MM-DD

## Standard Knockdown Pedal Exerciser For upper & lower body exercise

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