MoVeS

Wooden Incline Board Ladder:

REF

03-040001

UDI-DI

05420063003513

MVS In Motion
Westdijk 150
2830 Tiselt (Willebroek)
Belgium - Europe
www.mvs-in-motion.com
Find us on  & 

UK
CA

UK REP

CE MD

Qarad UK Ltd.
8 Northumberland Ave
Westminster, London WC2N 5BY
United Kingdom

ASYV3

en - Medical Device - Instructions for Use

Intended purpose: To tone, shape, stretch or rehabilitate muscle groups in lower legs or feet. **Intended user:** Lay person.

Indications: Medical conditions resulting in loss of strength, flexibility and coordination of lower legs, feet and ankle muscles. **Contra-indications:** In case of unusual swelling or discomfort, use should be discontinued.

Precautions: Consult a healthcare professional before use. Use on a stable, dry surface. Use of rubber soled shoe is recommended. Warm-up muscles before using device. Max user weight 150kg.

When assembling, always ensure that the tab on the incline platform is fully inserted into one of the four holes of the ladder section. Test the stability with one foot before stepping on with your full weight.

General information: Exercises should be done at a slow, comfortable pace.

Report any serious incident (death, serious deterioration of health, serious public health threat) to the Manufacturer and your Competent Authority.

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